ThermaD.A.S. D.A.S. ProPlasma & EndoD.A.S.

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Bibliografia

Long-wave plasma radiofrequency ablation for treatment of xanthelasma palpebrarum

Baroni A. – Journal of Cosmetic Dermatology 2018;001–3

Xanthelasma palpebrarum is the most common type of xanthoma affecting the eyelids. It is characterized by asymptomatic soft yellowish macules, papules, or plaques over the upper and lower eyelids.

This study aims to evaluate the effectiveness of RF ablation in the treatment of xanthelasma palpebrarum, with D.A.S. Medical portable device (Technolux, Italia), a radiofrequency tool working with long-wave plasma energy and without anesthesia.

Twenty patients, 15 female and 5 male, affected by xanthelasma palpebrarum, were enrolled for long-wave plasma radiofrequency ablation treatment. The treatment consisted of 3/4 sessions that were carried out at intervals of 30 days. Treatments were well tolerated by all patients with no adverse effects and optimal aesthetic results. The procedure is very fast and can be performed without anesthesia of the low because and tolerable pain stimulation. Long-wave plasma radiofrequency ablation is an effective option for treatment of xanthelasma palpebrarum and adds an additional tool to the increasing list of medical devices for aesthetic treatments.

Total Rejuvenation of the Face with Plasma Device

By Dr. Lebbar

PMFA News Vol.4, Ed.4 April-May 2017

The use of plasma energy is a new concept in anti-aging and aesthetic medicine. Since it has been introduced in the aesthetic field, many studies have been conducted to try and understand the effects on collagen synthesis which may explain the long lasting results.

The purpose of this study was to de safely with plasma energy and the ideal power settings of the plasma device to be used to have maximum stimulation of collagen synthesis.

The plasma dermabrasion technique is a very safe technique which can be used for the rejuvenation of the face: treating the eyelids, crow's feet, smoker's lines and it is also very capable as an electropeel.

The power and frequency settings of the medical device must be low while treating the eyelids to avoid the hyperpigmentation

secondary to excessive thermal damage. The perioral smoker's lines need a higher setting of power and frequency to obtain the maximum stimulation of collagen synthesis.

D.A.S. medical: the new plasma energy modulated in frequency By Dr. Lebbar

Ambulatorio Medico, nr. 52 Sept-Dec 2017

D.A.S. *medical* allows modulation of power and frequency. This means that the operator can decide, according to skin type, type of skin (thin or thick) and treatment, both quantity of tissue to remove and the thermal damage in a demined time on a specific area. Frequency modulation allows doctor to have total control on the treatment and to act in total safety being sure to avoid risks of secondary hyper-pigmentation on all skin types. D.A.S. medical device allows to perform different dermatological and medical aesthetic treatments, e.g. removal of dermal fibromas and eyelid skin laxity, in total safety. I have treated darker skin types without any risk of secondary hyper-pigmentation, since the device allows the modulation of frequency.

Total Rejuvenation of the Face with Plasma Device By Dr. Lebbar - AMWC 2017

TOTAL REJUVENATION OF THE FACE WITH PLASMA DEVICE

Dr Noura Lebbar, Cosmetic Physician, Milan, Italy

(mail dott.noura.lebbar@gmail.com)

Background / Introduction

The plasma energy is a new concept in Anti-Aging and Aesthetic Medicine. Since il has been introduced in the Aesthetic field, many studies have been conduced trying to understand the effects on the collagen synthesis which could explain the long lasting results.

The purpose of the study is to define the face areas to be treated safety with the plasma energy and the ideal power settings of the plasma device to be used to have maximum stimulation of collagen synthesis.

Material & Method

The medical device used is D.A.S. medical, is a plasma device which permits a controlled ionisation of the gas contained in the air producing plasma. This mechanism allows when the device is very closed to the skin, the sublimation of the upper layers of the epidermis. The stimulation of a new collagen synthesis is due to the controlled skin thermic damage (controlled thermolysis) due to the plasma energy. The medical plasma device has 5 levels of power settings and 5 levels of frequency setting. The level of frequency defines the depth of energy penetration in the skin, while the level of the power defines the degree of the thermic damage.

D.A.S. medical is the only plasma device with the possibility to settle the frequency in order to control the thermic damage and the risk of the hyperpigmentation.

I have treated 20 patients: 5 of them for the eyelids rejuvenation, 5 for the smoker lines, 5 for the crow's feet and 5 patients for the enlarged skin pores using plasma electropeel.

I have excluded tanned patients and patients with phototype 5 and more in the scale of fitzpatrick to avoid the risk of hyperpigmentation. No medical preparation was performed except an anesthetic cream 30.

Treatment of the upper & lower eyelids After 3 Months & 2 sessions

Perioral Area: After 3 weeks & only 1 session

No medical preparation was performed except an anesthetic cream 30 minutes before the procedure (lidocaine 30%). The session of plasma electropeel was followed always by an application of a mask with retinoid acid and growth factors. I have been operating using a triangle technique, i pulled small spots on the superficial skin layers, leaving 1 millimeter between every spot (healthy skin area), this method allows a faster recovery time. To treat the crow's feet, i have treated only the elevated skin from one part and the other of the wrinkle.

Results

75% of the patients were satisfied after only one session. I have performed 2 sessions for 2 patients treated for the crow's feet, 3 electropeel sessions for 4 patients, 2 sessions of eyelids rejuvenation for 2 patients, 3 sessions for smoker lines for 2 patients.

4 patients treated for the eyelids rejuvenation had an important skin tightening after only one session. The side effects were eyelids edema for 3 patients treated for the eyelids rejuvenation which lasts 5 days, 2 patients treated for the smoker lines had a redness lasting 7 days, 1 patient treated for the crow's feet had also an erythema lasting 10 days after the removal of the crust, 5 patients had a redness after the electropeel lasting 2 days. Only one patient treated for the eyelid rejuvenation had a hyper pigmentation secondary to an important eyelid inflammation not treated correctly.

Conclusion

The plasma dermabrasion technique is a very safe technique which can be used for the rejuvenation of the face: treating the eyelids, the crow's feet, the smoker lines and it is also very performant as an electropeel. The power and frequency settings of the medical device must be low while treating the eyelids to avoid the hyper pigmentation secondary to an excessive thermic damage. The smoker lines need a higher setting of power and frequency to obtain the maximum stimulation of collagen synthesis.

Treatment of the lower eyelids: 3 weeks after only 1 session

After 1 Month & only 1 session

La Dermo Ablation Surgery

by Dr. lozzo & Dr. Antonucci

Ambulatorio Medico, n.43 Sett-Dic 2014

The ever growing demand for minimally invasive treatments with rapid recovery has driven my interest towards a new technology called **Dermo Ablation Surgery (D.A.S.)** that has proved to be useful in the treatment of benign cutaneous lesions and of light or moderate skin eyelid laxity. The Dermo Ablation Surgery is an ablative long-waves radiofrequency that exploits a well-known principle, that of radiofrequency, that has been improved and enriched. D.A.S. technology is contained within a device consisting of a handpiece with a sterile metal electrode and a controller that allows the regulation of the working values. Approaching the electrode near the skin, it creates a voltaic arch that generates an electromagnetic wave responsible for the effect on the skin itself.

In view of an improvement and rejuvenation of the face, it is possible to associate the D.A.S. *medical* to other treatments such as botulinum toxin or hyaluronic acid.